


# Complete catalogue

Home care

Perfectly suited for care.

.bock<sup>®</sup>///



“Healthcare is one of the most difficult arts. It is motivated by compassion but driven by knowledge alone.”

Mary Adelaide Nutting, first professor of nursing, 1925

A health bed from bock is always the outcome of synergy. Our internal research and development centre forms the basis for this. Here skills from the engineering, scientific, medical and healthcare domains are combined. A care bed needs to satisfy many requirements and personal needs.

So whenever a new development begins, we aim to achieve an ideal mix of individuality, functionality and economy.

After all, our aspiration can only culminate in visible results in everyday care – satisfaction on the part of operators and users alike – when all these issues have been taken into account.

Yours


A handwritten signature in black ink, reading 'Klaus Bock'. The script is fluid and cursive, with the first letters of 'Klaus' and 'Bock' being capitalized and prominent.

Klaus Bock



# Contents

<b>Bock as an innovator</b>	<b>4</b>	<b>Technology and safety</b>	<b>36</b>
<b>Quality made in Germany</b>	<b>5</b>	<b>ripolux® neo</b>	<b>37</b>
<b>Home care</b>	<b>6</b>	<b>Accessories</b>	<b>38</b>
<b>Special beds</b>	<b>7</b>	Accessories for more individuality	39
<b>Overview</b>	<b>8</b>	<b>Bed side cabinets</b>	<b>40</b>
domiflex 2 with Push-and-Ready	10	<b>Option management</b>	<b>42</b>
domiflex	12	Lying surfaces	42
belluno	14	Side rails	43
combiflex nursing care bed	16	<b>Mattresses</b>	<b>44</b>
belluno bibs (bed-in-bed system)	18	<b>New releases 2016</b>	<b>46</b>
combiflex bibs (bed-in-bed system)	20	dino   Cot & adolescent bed	46
domiflex classic low   Low bed	22	domiflex 2 wash with Push-and-Ready	48
domiflex low   Low bed	24	<b>Staying, living, caring.</b>	<b>50</b>
practico ultra low 9.5/80   Ultra low bed	26		
domiflex 185	28		
adiflex 220   Bariatric bed	30		
adilec 220   Bariatric bed	32		
adilec 280   Bariatric bed	34		



## bock as an innovator

"For us at Hermann Bock GmbH, technical progress does not count as innovation until it helps to improve the quality of life of people in need of healthcare, and consequently increases their zest for life. That is how we think and act when developing high-grade health beds. We use our long-standing experience and technical know-how to create products which are geared to the individual needs of people in care and the requirements of nursing personnel, and which extend well beyond the standards set out by the law. Due to our working processes beds with added values are the result, extra convenience and comfort and a greater sense of well-being under economical aspects."







## Quality– Made in Germany

More than ninety years of successful company history, more than ninety years of unstoppable innovation capacity characterize Hermann Bock GmbH. Numerous forward-thinking developments and patented solutions have put us in a visionary position in the healthcare market. Our special interest in lying healthy began with the production of slatted frames in the early 1970s. Development and manufacture of a dual-motor slatted frame was followed by the first electrically height-adjustable care bed in 1983. Since 1992, Hermann Bock has offered special furnishing concepts for nursing and healthcare institution and meanwhile became one of the leading German manufacturers of beds for therapy and care. At our plant in Verl approximately 120 employees are working continuously on the enhancement and improvement of the health care products, from preventative lying surface systems to complete furnishing concepts.

With ideas, knowledge and modern production techniques we realize our high requirements on quality, optimal safety and comfort.

With every innovation we want to contribute to more quality of life by increased well-being and an efficient health care system. This is our aim now and in future.

High quality and economical products go hand in hand at bock.

At our own innovation centre, our engineers develop bed systems that are in the line with market and personal requirements.

Prototypes are tested and test methods are established to analyze the risks. Modern machinery for punching, bending, robotic welding and electrostatic powder-coating are available for serial production at our manufacturing facilities that extend over an area of 9000 square meters.

# Home care

In their unique quality, Hermann Bock therapy and rehabilitation beds unite the wealth of knowledge that has been acquired from years of experience as a manufacturer with the latest research solutions, science and practice. Our healthcare products equally meet the patient's needs for maximum comfort and dignified living.

The highest safety standards and first-class functionality have proved themselves in professional healthcare, day after day, since years.

Our future efforts in co-operation with experts in rehabilitation and furniture sectors will help to develop new healthcare products and enhance existing ones, and thus continue to provide optimum support for quality of life.





## Special beds

Special needs require special measures. Our long experience in dealing with physical disabilities and restriction as a consequence of illness and accidents have repeatedly confirmed that, it is impossible to generalize people who are in need of therapy. Therefore an important aim of our ambitions is to develop special bed systems, corresponding to the requirements of people, for whom; due to their physical condition or diseases,

the standard care system does not provide the optimal care and treatment. With our special beds, we have laid the foundation to provide the support for the recuperation process while creating professional care conditions at the same time.

# Overview

## Our standard beds



► **domiflex 2 with Push-and-Ready**  
– Page 10



► **domiflex**  
– Page 12



► **belluno**  
– Page 14



► **combiflex**  
– Page 16



► **belluno bibs (bed-in-bed system)**  
– Page 18



► **combiflex bibs (bed-in-bed system)**  
– Page 20



## Our special beds



- ▶ **domiflex classic low**
- Low bed
- Page 22



- ▶ **domiflex low**
- Low bed
- Page 24



- ▶ **practico ultra low 9.5/80**
- Ultra low bed for a person weight up to 185 kg
- Page 26



- ▶ **domiflex 185**
- Bariatric bed for a person weight up to 185 kg
- Page 28



- ▶ **adiflex 220 | Bariatric bed**
- Bariatric bed for a person weight up to 220 kg
- Page 30



- ▶ **adilec 220 | Bariatric**
- Bariatric bed for a person weight up to 220 kg
- Page 32

## New releases 2016



- ▶ **adilec 280 | Bariatric bed**
- Bariatric bed for a person weight up to 280 kg
- Page 34



- ▶ **dino cot & adolescent nursing care bed**
- Page 46



- ▶ **domiflex 2 wash with Push-and-Ready**
- Page 48



# domiflex 2 with Push and Ready

Perfectly equipped for home care, with the Push-and-Ready connection technique



## Specifications – domiflex 2

Total weight:	approx. 83 kg
Lying surface dimension:	90 x 200 cm
External dimension:	103 x 213 cm
Safe working load:	190 kg
Maximum person weight:	155 kg
Height adjustment:	39 – 80 cm
Lifter space:	> 15 cm
Side rail height:	37 cm with wooden slats 34 cm with ripolux® neo
Special dimensions:	220 cm with bed extension only In a special width only available in: 100 cm as domiflex 185 120 cm as adiflex 220

## Motors

Water-resistant to IPX4, high quality spiral cable with pull relief and anti-kink device, as well as a safety extra-low voltage system.

## Features

### Art.-No.

271	Wooden side rails	90 x 200 cm
271	Steel side rails	90 x 200 cm

Subject to technical modification.



### 4-section

Electrical adjustment of the lying surface including the automatic triple knee break action.



3 Beds on a pallet of 120x80 cm

## Advantages at a glance

- Assembly in less than 5 minutes without any tools (by only one technician)
- Time saving innovative Push-and-Ready™ connection technique
- No loose parts required for assembly due to pre-installed components
- Simple handling thanks to the light weight of all components
- Space-saving storage and effective transportation (three beds on one pallet)
- Completely integrated van-transport system
- Easy cleaning thanks to smooth surfaces
- Increased safe capacity of 190 kg and maximum person weight of 155 kg
- Guarantee on motors extended to 5 years

### Basic equipment:

Safety extra-low voltage system  
Lifting pole with hand grip  
Colour: optionally in D7935 Havana or F37/009 Steamed beech  
Wooden slats  
Integrated steel or wooden side rails  
Push-and-Ready connection technique  
Transport and storage system  
Integrated Van-transport system (Only with the universal transport system)

### Extras:



Bed extension for the length to 220 cm  
Covering of the motor  
Comfort-sitting position/ Reverse Trendelenburg  
Mattress compensation  
ripolux® neo-spring lying surface system incl. mattress compensation  
Steel lying surface with removable bars  
Full covering of the end panels in D7935 Havana or F37/009 Steamed beech

### Accessories:

Infusion holder  
Mattresses  
Mattress extension for the length to 220 cm  
Side rail attachment  
Side rail bumpers  
Side rail bumpers for side rails with attachment  
Light underneath the bed  
And many more

### Safety guarantees:

- Primary fuse in the power plug
- Secondary fuse/polyswitch
- Thermal release
- Locking function on the hand control
- Hand control with single faulty safety
- Special spiral mains cable
- Pull relief at power plug and anti-kink device
- Water-resistant to IPX4
- Distortion-free, highly stable side rails





# domiflex

The well-known classic

## Specifications – domiflex

Total weight:	approx. 82 kg
Lying surface dimension:	90 x 200 cm
External dimension:	103 x 221 cm
Safe working load:	170 kg
Maximum person weight:	135 kg
Height adjustment:	40 – 81 cm
Lifter space:	> 15 cm
Side rail height:	37 cm with wooden slats 34 cm with ripolux® neo
Special dimensions:	180 – 220 cm length In a special width only available in: 100 cm as domiflex 185 120 cm as adiflex 220

## Motors

Water-resistant to IPX4, high quality spiral cable with pull relief and anti-kink device, as well as a safety extra-low voltage system.

## Features

### Art.-No.

261	Wooden side rails	90 x 200 cm
261	Steel side rails	90 x 200 cm

Subject to technical modification.



### 4-section

Electrical adjustment of the lying surface including the automatic triple knee break action.



Transport system for domiflex nursing care bed



### Basic equipment:

Safety extra-low voltage system  
Lifting pole with hand grip  
Colour: optionally in D7935 Havana or F37/009 Steamed beech  
Wooden slats  
Integrated steel or wooden side rails  
Transport and storage system

### Extras:


Bed extension for the length to 220 cm  
Covering of the motor (only available in Havana colour)  
Comfort-sitting position/ Reverse Trendelenburg  
Mattress compensation  
ripolux® neo- spring lying surface system incl. mattress compensation  
Steel lying surface with removable bars  
Special length

### Accessories:

Bed lamps  
Infusion holder  
Mattresses  
Mattress extension for the length to 220 cm  
Side rail bumpers  
Side rail bumpers for side rails with attachment  
Side rail attachment  
And many more

### Safety guarantees:

- Primary fuse in the power plug
- Secondary fuse/polyswitch
- Thermal release
- Locking function on the hand control
- Hand control with single faulty safety
- Special spiral mains cable
- Pull relief at power plug and anti-kink device
- Water-resistant to IPX4
- Distortion-free, highly stable side rails





# belluno

A modern classic for proper rest

## Specifications – belluno

Total weight:	approx. 98 kg
Lying surface dimension:	90 x 200 cm
External dimension:	103 x 213 cm
Safe working load:	200 kg
Maximum person weight:	165 kg
Height adjustment:	36 – 80 cm
	29 – 73 cm as low version
Lifter space:	> 15 cm
Side rail height :	37 cm with wooden slats
	34 cm with ripolux® neo
Special dimensions:	80 – 140 cm width
	180 – 220 cm length
	Low version not available in
	80 or 110 cm lying surface width

## Motors

Water-resistant to IPX4, high quality spiral cable with pull relief and anti-kink device, as well as a safety extra-low voltage system.

## Features

### Art.-No.

281	Wooden side rails	90 x 200 cm
281	Steel side rails	90 x 200 cm

Subject to technical modification.



### 4-section

Electrical adjustment of the lying surface including the automatic triple knee break action.



Transport system for belluno nursing care bed

### Basic equipment:

Safety extra-low voltage system  
Colour: optionally in Havana D7935 or steamed Beech F37/009  
Selected special colours are possible against extra charge<sup>2</sup>  
Lifting pole with hand grip  
Wooden slats  
Integrated steel or wooden side rails

### Extras:

Leg position parallel to lying surface with 90° end position  
Bed extension for the length to 220 cm  
Wooden slats with tension regulation  
Mattress compensation  
Integrated head rest for 5-section lying surface<sup>1</sup>  
Low version  
ripolux® neo- spring lying surface system incl. mattress compensation  
Steel lying surface with removable bars  
Special width/ Special length  
Central braking

### Accessories:

Bed lamps  
Infusion holder  
Mattresses  
Mattress extension for the length to 220 cm  
Side rail bumpers  
Side rail bumpers for side rails with attachment  
Side rail attachment  
And many more

### Safety guarantees:

- Primary fuse in the power plug
- Secondary fuse/polyswitch
- Thermal release
- Locking function on the hand control
- Hand control with single faulty safety
- Special spiral mains cable
- Pull relief at power plug and anti-kink device
- Water-resistant to IPX4
- Distortion-free, highly stable side rails

<sup>1</sup> Inclusive mattress compensation and wooden slats with tension regulation

<sup>2</sup> The special colours can be found on page 51







# combiflex

The changeable nursing care bed, 220 kg safe working load

## Specifications – combiflex nursing care bed

Total weight:	approx. 113 kg
Lying surface dimension:	90 x 200 cm
External dimension:	103 x 213 cm
Safe working load:	220 kg
Maximum person weight:	185 kg
Height adjustment:	37 – 81 cm
Lifter space:	> 15 cm
Side rail height :	37 cm with wooden slats 34 cm ripolux® neo
Special dimensions:	80 – 140 cm width 190 – 220 cm length


## Motors

Water-resistant to IPX4, high quality spiral cable with pull relief and anti-kink device, as well as a safety extra-low voltage system.

## Features

293	Wooden side rails	90 x 200 cm
293	Steel side rails	90 x 200 cm

Subject to technical modification.



## 4-section

Electrical adjustment of the lying surface including the automatic triple knee break action.



Transport system for combiflex nursing care bed





### Basic equipment:

Safety extra-low voltage system  
775 mm head bar  
Lifting pole with hand grip  
Colour: optionally in Havana D7935 or steamed Beech F37/009  
Selected special colours are possible against extra charge<sup>1</sup>  
Wooden slats with tension regulation  
Wooden end panel type VIII  
Integrated steel or wooden side rails  
Comfort-sitting position/ Reverse Trendelenburg  
Mattress compensation  
Hand control with locking function

### Extras:

Leg position parallel to lying surface with 90° end position  
Integrated head rest for 5-section lying surface  
ripolux® neo- spring lying surface system incl. mattress compensation  
Steel lying surface with removable bars  
Special width/ Special length  
Steel lying surface with removable bars  
Central braking

### Accessories:

Bed lamps  
Bed extension for the length to 220 cm  
Infusion holder  
Mattresses  
Mattress extension for the length to 220 cm  
Side rail bumpers  
Side rail bumpers for side rails with attachment  
Side rail attachment  
Transport and storage system  
And many more

### Safety guarantees:

- Primary fuse in the power plug
- Secondary fuse/polyswitch
- Thermal release
- Locking function on the hand control
- Hand control with single faulty safety
- Special spiral mains cable
- Pull relief at power plug and anti-kink device
- Water-resistant to IPX4
- Distortion-free, highly stable side rails

<sup>1</sup> The special colours can be found on page 51





# belluno bibs (bed-in-bed system)

With bed insert everywhere at home

## Specifications – belluno bibs (Bed-in-bed system)

Total weight:	approx. 65 kg
Lying surface dimension:	90/100 <sup>1</sup> x 190 or 200 cm
Nominal dimension:	89 x 187/197 cm
Safe working load:	200 kg
Maximum person weight:	165 kg
Height adjustment:	31 – 75 cm
Lifter space:	11 cm
Special dimensions:	80 – 140 cm width 180 – 220 cm length

## Motors


Water-resistant to IPX4, high quality spiral cable with pull relief and anti-kink device, as well as a safety extra-low voltage system.

## Features


### Art.-No.

282	90/100 x 190 cm
282	90/100 x 200 cm

Subject to technical modification.



Base in a length of 1240 mm  
(Only upon request)



Base in a length of 1560 mm,  
standard base width

## 4-section

Electrical adjustment of the lying surface including the automatic triple knee break action.

Transport system for belluno bibs

### Basic equipment:

Safety extra-low voltage system  
Wooden slats  
Provided feet with balance screw for free standing assembly  
Mattress brackets for width extension

### Extras:

Leg position parallel to lying surface with 90° end position  
Wooden slats with tension regulation  
Higher feet to increase the height by approx. 4 cm  
Mattress compensation  
Integrated head rest for 5-section lying surface  
ripolux® neo- spring lying surface system incl. mattress compensation  
Steel lying surface with removable bars  
Mattress bracket on the foot end


### Accessories:

Lifting pole with hand grip  
Companion bed 90/100 cm width (wood work only)  
Bed lamps  
Foot bar 190 cm for 90 cm width (to extend to 200 cm)  
Foot bar 200 cm for 90 cm width (to shorten to 190 cm)  
Infusion holder  
Mattresses  
Side rail bumpers  
Side rails for bed insert, 100 cm fitting  
Side rails for bed insert, 90 cm fitting  
Transport and storage system

### Safety guarantees:

- Primary fuse in the power plug
- Secondary fuse/polyswitch
- Thermal release
- Locking function on the hand control
- Hand control with single faulty safety
- Special spiral mains cable
- Pull relief at power plug and anti-kink device
- Water-resistant to IPX4

<sup>1</sup> With the aid of the mattress bracket extensible in width.





# combiflex bibs (Bed-in-bed system)

The modular system, 220 kg safe capacity

## Specifications – combiflex bibs (Bed-in-bed system)

Total weight:	approx. 85 kg
Lying surface dimension:	90/100 <sup>1</sup> x 190 or 200 cm
Nominal dimension:	89 x 187/197 cm
Safe working load:	220 kg
Maximum person weight:	185 kg
Height adjustment:	33 – 77 cm
Lifter space:	11 cm
Special dimensions:	80 – 140 cm width 190 – 220 cm length


## Motors

Water-resistant to IPX4, high quality spiral cable with pull relief and anti-kink device, as well as a safety extra-low voltage system.

## Modular system

<b>Art.-No. 291</b>	combiflex bed insert
<b>Art.-No. 292</b>	Wooden bed frame (end panels, side rails, lifting pole, castors)
<b>Art.-No. 293</b>	combiflex nursing care bed (bed insert, wooden bed frame, lifting pole, castors)
<b>Art.-No. 90753</b>	Transport and storage system

Subject to technical modification.




## 4-section

Electrical adjustment of the lying surface including the automatic triple knee break action.



Transport system for combiflex bibs





### Basic equipment:

Safety extra-low voltage system  
775 mm head bar  
Wooden slats with tension regulation  
Comfort-sitting position / Reverse Trendelenburg  
Mattress compensation  
Mattress brackets for width extension  
Hand control with locking function  
Feet incl. balancing screws

### Extras:

Leg position parallel to lying surface with 90° end position  
Feet increase by approx. 4 cm  
Mattress bracket on the foot end  
Integrated head rest for 5-section lying surface  
ripolux® neo- spring lying surface system incl. mattress compensation  
Steel lying surface with removable bars  
Special length / Special width


### Accessories:

Lifting pole with hand grip  
Bed lamps  
Foot bar 190 cm for 90 cm width (to extend to 200 cm)  
Foot bar 200 cm for 90 cm width (to shorten to 190 cm)  
Wooden end panel type VIII  
Infusion holder  
Mattresses  
Mattress extension for the length to 220 cm  
Side rails for bed insert, 100 cm fitting  
Side rails for bed insert, 90 cm fitting  
Transport and storage system

### Safety guarantees:

- Primary fuse in the power plug
- Secondary fuse/polyswitch
- Thermal release
- Locking function on the hand control
- Hand control with single faulty safety
- Special spiral mains cable
- Pull relief at power plug and anti-kink device
- Water-resistant to IPX4

<sup>1</sup> With the aid of the mattress bracket extenable in the width.





# domiflex classic low

Low entrance height, high quality

## Specifications – domiflex classic low

Total weight:	approx. 82 kg
Lying surface dimension:	90 x 200 cm
External dimension:	103 x 221 cm
Safe working load:	170 kg
Maximum person weight:	135 kg
Height adjustment:	21 – 61 cm
Lifter space:	> 15 cm
Side rail height :	37 cm with wooden slats 34 cm with ripolux® neo
Special dimensions:	180 – 220 cm length Special widths not available

## Motors

Water-resistant to IPX4, high quality spiral cable with pull relief and anti-kink device, as well as a safety extra-low voltage system.

## Features

### Art.-No.

263	Wooden side rails	90 x 200 cm
263	Steel side rails	90 x 200 cm

Subject to technical modification.



### 4-section

Electrical adjustment of the lying surface including the automatic triple knee break action.



Transport system for domiflex classic low

### Basic equipment:

Safety extra-low voltage system  
Colour: D 7935 Havana  
Motor covering of the actuators in matching colour  
Wooden slats  
Integrated steel or wooden side rails  
Transport and storage system

### Extras:

Bed extension for the length to 220 cm  
Trendelenburg / Reverse Trendelenburg  
Mattress compensation  
ripolux® neo- spring lying surface system incl. mattress compensation  
Steel lying surface with removable bars  
Special length

### Accessories:

Lifting pole with hand grip  
Bed lamps  
Infusion holder  
Mattresses  
Mattress extension for the length to 220 cm  
Side rail bumpers  
Side rail bumpers for side rails with attachment  
Side rail attachment  
And many more

### Safety guarantees:

- Primary fuse in the power plug
- Secondary fuse/polyswitch
- Thermal release
- Locking function on the hand control
- Hand control with single faulty safety
- Special spiral mains cable
- Pull relief at power plug and anti-kink device
- Water-resistant to IPX4
- Distortion-free, highly stable side rails





# domiflex low

Very low - highly protective

## Specifications – domiflex low

Total weight:	approx. 80 kg
Lying surface dimension:	90 x 200 cm
External dimension:	103.5 x 215.5 cm
Safe working load:	170 kg
Maximum person weight:	135 kg
Height adjustment:	24 – 64 cm
Lifter space:	> 15 cm
Side rail height :	37 cm with wooden slats 34 cm ripolux® neo
Special dimensions:	180 – 220 cm length Special widths not available

## Motors

Water-resistant to IPX4, high quality spiral cable with pull relief and anti-kink device, as well as a safety extra-low voltage system.

## Features

### Art.-No.

266

Wooden side rails 90 x 200 cm

Subject to technical modification.



### 4-section

Electrical adjustment of the lying surface including the automatic triple knee break action.



### Basic equipment:

Safety extra-low voltage system  
Colour: optionally Havana D7935 or steamed Beech F37/009  
Wooden slats  
Integrated steel or wooden side rails  
Low entrance height  
Full covering of the end panels

### Extras:


Bed extension for the length to 220 cm  
Hand control with Reverse-Trendelenburg  
Mattress compensation  
ripolux® neo- spring lying surface system incl. mattress compensation  
Steel lying surface with removable bars

### Accessories:

Roll-down mattress for low bed  
Lifting pole with hand grip  
Bed lamps  
Infusion holder  
Mattresses  
Mattress extension for the length to 220 cm  
Side rail bumpers  
Side rail bumpers for side rails with attachment  
Side rail attachment  
Transport and storage system

### Safety guarantees:

- Primary fuse in the power plug
- Secondary fuse/polyswitch
- Thermal release
- Locking function on the hand control
- Hand control with single faulty safety
- Special spiral mains cable
- Pull relief at power plug and anti-kink device
- Water-resistant to IPX4
- Distortion-free, highly stable side rails





# practico ultra low 9.5/80

The barrier-free bed - for independent mobility and autonomy

## Specifications – practico low 9.5/80

Lying surface dimension:	90 x 200 cm
External dimension:	101.5 x 240 cm
Safe working load:	220 kg
Maximum person weight:	185 kg
Height adjustment:	9.5 – 80 cm
Lifter space:	> 15 cm
Side rail height :	39.5 cm with wooden slats, 35 cm with ripolux® neo


## Motors

Water-resistant to IPX4, high quality spiral cable with pull relief and anti-kink device, as well as a low voltage system.

## Features

Art.-No.	Lying surface dimension	Height adjustment
422	90 x 200 cm	9.5 to 80 cm

Subject to technical modification.



## 4-section

Electrical adjustment of the lying surface including the automatic triple knee break action.

### Basic equipment:

Safety extra-low voltage system  
Colour: F37/009 Steamed beech  
Type 32 end panel  
Wooden slats with tension regulation  
Integrated bed extension to 220 cm  
Comfort-sitting position / Trendelenburg  
Mattress compensation  
Ultra low entrance height  
Lighting underneath the lying surface  
Castors with single brakes and integrated wall protection system

### Extras:

End panel type 31  
End panel type 33  
Bed extension for the length to 220 cm  
Aluminium lying surface with removable bars  
ripolux® neo- spring lying surface system incl. mattress compensation  
Steel lying surface with removable bars  
Continuous side rails  
Long divided side rails (bed can only be lowered to 15.5 cm)

### Accessories:

Lifting pole with hand grip  
Bed lamps  
Mattresses  
Mattress extension for the length to 220 cm  
Side rail bumpers  
And many more

### Safety guarantees:


- Primary fuse in the power plug
- Thermal release
- Locking function on the hand control
- Hand control with single faulty safety
- Special spiral mains cable
- Pull relief at power plug and anti-kink device
- Water-resistant to IPX4
- Distortion-free, highly stable side rails



End panel type 32



End panel type 31







# domiflex 185

The low-cost variant for a person weight up to 185 kg

## Specifications – domiflex 185 kg

Total weight:	approx. 93 kg
Lying surface dimension:	90 x 200 cm
External dimension:	103 x 214 cm
Safe working load:	220 kg
Maximum person weight:	185 kg
Height adjustment:	40 – 81 cm
Lifter space:	> 15 cm
Side rail height :	37 cm with wooden slats 34 cm with ripolux® neo
Special dimensions:	100 cm width 180 – 220 cm length 120 cm as adiflex 220

## Motors

Water-resistant to IPX4, high quality spiral cable with pull relief and anti-kink device, as well as a safety extra-low voltage system.

## Features

### Art.-No.

264	Wooden side rails	90/100 x 200 cm
264	Steel side rails	90/100 x 200 cm

Subject to technical modification.



### 4-section

Electrical adjustment of the lying surface including the automatic triple knee break action.



Transport system for domiflex 185



### Basic equipment:

Safety extra-low voltage system  
Lifting pole with hand grip  
Motor covering of the actuators in matching colour  
Colour: Havana D 7935  
Wooden slats with tension regulation  
Integrated steel or wooden side rails  
Transport and storage system

### Extras:




Bed extension for the length to 220 cm  
Mattress compensation  
ripolux® neo- spring lying surface system incl. mattress compensation  
Steel lying surface with removable bars  
Special length

### Accessories:

Bed lamps  
Infusion holder  
Mattresses  
Mattress extension for the length to 220 cm  
Side rail bumpers  
Side rail bumpers for side rails with attachment  
Side rail attachment  
And many more

### Safety guarantees:

- Primary fuse in the power plug
- Secondary fuse/polyswitch
- Thermal release
- Locking function on the hand control
- Hand control with single faulty safety
- Special spiral mains cable
- Pull relief at power plug and anti-kink device
- Water-resistant to IPX4
- Distortion-free, highly stable side rails





# adiflex 220

Comfort for heavy patients up to 220 kg

## Specifications – adiflex 220 Bariatric bed

Total weight:	approx. 126 kg
Lying surface dimension:	120 x 200 cm
External dimension:	133 x 217 cm
Safe working load:	255 kg
Maximum person weight:	220 kg
Height adjustment:	23 – 83 cm
Lifter space:	> 15 cm
Side rail height :	37 cm with wooden slats


## Motors

Water-resistant to IPX4, high quality spiral cable with pull relief and anti-kink device, as well as a safety extra-low voltage system.

## Features

Art.-No.	Lying surface dimension	External dimension
267	120 x 200 cm	133 x 217 cm
267	120 x 220 cm	133 x 237 cm

Subject to technical modification.



## 4-section

Electrical adjustment of the lying surface including the automatic triple knee break action.

### Basic equipment:



Safety extra-low voltage system  
Colour: Havana D 7935  
Wooden slats with tension regulation  
Integrated wooden side rails

### Accessories:

Roll-down mattress  
Lifting pole with hand grip  
Bed lamps  
Infusion holder  
Mattresses  
Side rail bumpers  
Side rail bumpers for side rails with attachment  
Special mattress

### Safety guarantees:

- Primary fuse in the power plug
- Secondary fuse/polyswitch
- Thermal release
- Locking function on the hand control
- Hand control with single faulty safety
- Special spiral mains cable
- Pull relief at power plug and anti-kink device
- Water-resistant to IPX4
- Distortion-free, highly stable side rails





# adilec 220

Premium nursing care bed for a person weight up to 220 kg

## Specifications – adilec 220 Bariatric bed

Total weight:	approx. 155 kg
Lying surface dimension:	look features
External dimension:	look features
Safe working load:	255 kg
Maximum person weight:	220 kg
Height adjustment:	40 – 81 cm
Lifter space:	> 15 cm
Side rail height :	37 cm

## Motors

Water-resistant to IPX4, high quality spiral cable with pull relief and anti-kink device as well as safety in power plug, mains isolation system and emergency lowering

## Features

Art.-No.	Lying surface dimension	External dimension
246	100 x 200 cm	116 x 223 cm
246	100 x 220 cm	116 x 243 cm

Subject to technical modification.



## 4-section

Electrical adjustment of the lying surface including the automatic triple knee break action.



### Basic equipment:



Colour: F37/009 Steamed beech  
Wooden slats with tension regulation  
Integrated wooden side rails  
Full covering of the end panels

### Accessories:

Lifting pole with hand grip  
Bed lamps  
Infusion holder  
Mattresses  
Side rail bumpers  
Side rail bumpers for side rails with attachment  
Special mattresses  
Central braking 2+2 system  
And many more

### Safety guarantees:

- Mains isolation
- Primary fuse in the power plug
- Secondary fuse/polyswitch
- Thermal release
- Locking function on the hand control
- Hand control with single faulty safety
- Special spiral mains cable
- Pull relief at power plug and anti-kink device
- Water-resistant to IPX4
- Distortion-free, highly stable side rails





# adilec 280

Robust premium care bed for a person weight up to 280 kg

## Specifications – adilec 280 Bariatric bed

Total weight:	approx. 210 kg
Lying surface dimension:	look features
External dimension:	look features
Safe working load:	315 kg
Maximum person weight:	280 kg
Height adjustment:	40 – 81 cm
Lifter space:	> 15 cm
Side rail height :	39 cm with wooden slats

## Motors

Water-resistant to IPX4, high quality spiral cable with pull relief and anti-kink device as well as safety in power plug, mains isolation system and emergency lowering.

## Features

Art.-No.	Lying surface dimension	External dimension
245	100 x 200 cm	116 x 223 cm
245	100 x 220 cm	116 x 243 cm
245	120 x 200 cm	136 x 223 cm
245	120 x 220 cm	136 x 243 cm
245	140 x 200 cm	156 x 223 cm
245	140 x 220 cm	156 x 243 cm

Subject to technical modification.



## 4-section

Electrical adjustment of the lying surface including the automatic triple knee break action.

### Basic equipment:


Colour: F37/009 Steamed beech  
Integrated wooden side rails  
Full covering of the end panels  
Central braking 2+2 system

### Accessories:

Lifting pole with hand grip  
Bed lamps  
Infusion holder  
Mattresses  
Side rail bumpers  
Side rail bumpers for side rails with attachment  
Special mattresses

### Safety guarantees:

- Mains isolation
- Primary fuse in the power plug
- Secondary fuse/polyswitch
- Thermal release
- Locking function on the hand control
- Hand control with single faulty safety
- Special spiral mains cable
- Pull relief at power plug and anti-kink device
- Water-resistant to IPX4
- Distortion-free, highly stable side rails





# Technology and safety

## **Technology and safety**

Every single bed that we manufacture is subjected to comprehensive functional tests here at our production site; this is part of our quality assurance system. Whenever our beds are delivered to a customer for the very first time, the person in charge chosen by the operator is trained to handle the beds properly. In addition Hermann Bock offers at any time training for the proper use and safety of the bed.

## **Manufacturer, operator and user – a team**

Safety regulations and measures, the manufacturer, operator and user form a unit, in which the responsibilities are clearly defined. The operator is any individual or juristic person that uses Hermann Bock care beds or has used them on their behalf, including specialty retailers, health insurances, care homes or other institutions. They are responsible for teaching the users how the beds are properly used and maintained. Users are all persons who have the necessary training and experience to use electric health beds from Hermann Bock.

Users are for example, nursing staff and relatives who take care of a family member. The manufacturer provides the user with the user's manual for the personal information and transmission. These documents are for the operator's own information and for the transmission to the user.

## **Safety inspections with service**

Inspections are conducted to keep the very high safety level that we aspire to, and are therefore an important measure to protect individuals. Medical equipments must be regularly inspected for safety at the intervals stated by the manufacturer and generally accepted guidelines. In everyday practice, safety measures are subject to different requirements and demands and furthermore to potential signs of wear. To avoid dangers it is essential to carry out the safety inspection in the line with national guidelines. The manufacturer has no influence on the extent to which the operator follows the rules. By offering time-saving services, bock makes it easier for you to follow the necessary safety rules.

A bed is only released for series production when the remaining risks are minimal. All health beds from Hermann Bock can be categorized as active medical products belonging to class 1 (lowest risk class) according to the rules 1 and 12 (Appendix IX of EU directive 93/42/EWG). Conformity with EU directive 93/42/EWG is ensured by means of a conformity assessment process.

As we are a manufacturer of medical products, our quality management system based on DIN EN ISO 90001 and DIN EN ISO 13485:2012, is an effective instrument to optimize the process at our company. By fulfilling this standard, we can assure our customers that our products meet legal requirements and therefore increase the satisfaction level of operators and users.





On the outside lying  
spring bows with a  
central spring bridge.

Ideal ventilation thanks to  
an open support surface.

## ripolux® neo

One spring support system,  
perfect for everyone.

The flexible lying surface system can be adjusted to fit the body's pressure distribution perfectly – for everyone, no matter how large or small, heavy or light, impaired or physically able. ripolux® neo ensures comfortable, ergonomic positioning in bed, in every sitting and lying position, and for every medical need. When patients' senses and/or bodies are affected by disease or paralysis, ripolux® neo helps – with selective pressure relief plus targeted back and body support. With ripolux® neo, carers can relieve pressure on patients' backs in no time at all – and support other areas of the body at the same time.

### The 1st choice in rehabilitation and care – ripolux® neo for all needs

- High-performance plastic for durable elasticity and a long lifetime
- Springs in three different degrees of tension for individual configuration
- Spring configuration individually adoptable without any tools
- Harmonious pressure balance thanks to an independent spring deflection and a very detailed arrangement of springs
- No moisture absorption and easy cleaning – for maximum hygiene
- Ideal ventilation of the mattress and a minimum of support surfaces – ideal mould protection
- The mattress profits from freedom of movement: it does not thread, get stuck, tilt or cause wrinkles



# Accessories

for more individuality



► Attachable bed lamp, model F,  
switchable via hand control




► Bed lamp underneath the  
lying surface, model UB,  
with motion sensor

Alternative:  
Bed lamp underneath the  
lying surface, model UF,  
switchable by hand control



► Attachable bed lamp, model E



► Stand-alone lifting pole



► Lifting pole with hand grip



► Grab rail



► Infusion holder duplex



► Infusion holder, height adjustable



► Infusion holder for lifting pole



► Attachable tray  
– to be placed on the side rails




► Universal clamp  
with hand control holder



► Universal clamp  
with urine bottle holder



► Universal clamp  
with drainage bottle holder



► Roll-down mattress for low beds



► Side rail attachment for continuous side rails  
(the side rail height will be increased by approx. 10 cm)



► Side rail bumpers





► Transport system  
for belluno and combiflex



► Transport system  
for bed-in-bed system


# Bedside cabinets

[Available D 7935 Havana and F 37/009 Steamed Beach]





- Standard 1300
- With drawer and door

- Standard 1301
- With drawer and deep drawer



- Standard 1302
- With drawer and deep drawer
- Shelf board



- Small server
- Continuously adjustable
- Swiveling tray

- Server 300
- Continuously adjustable
- Fixed deposit
- Swiveling tray with stop bars






- domiflex
- With drawer and door
- Height adjustable fold-away tray
- Tray attachable on the right and left side



- perfekta
- Integrated fold-away tray
- Drawer usable from both sides
- Open compartment
- 1 shelf behind door
- Doors on both sides



- verona 310
- Integrated tray with supporting lift, height adjustable
- Drawer usable from both sides
- Open compartment
- Doors on both sides
- 1 shelf behind door




- aluna server 351
- Tray available in different colours, insertion profile only alu colour
- Oval C-profile tube with 75 mm double steering castors
- Tray system (non fold-away) with height adjustment aid
- With fixed deposit and swiveling compact tray incl. stop bars
- Tray with automatic release, height adjustable by pneumatic spring
- With multifunctional bar



- aluna server 350
- With fold-away tray system
- Only available in alu colour

# Option management


## Lying surfaces



- ▶ Wooden slats
- Standard equipment



- ▶ Lying surface
- With removable steel bars



- ▶ ripolux® neo
- Spring system based on removable aluminium bars

## Side rails



► Wooden side rails




► Steel side rails




► Side rails for bed-in-bed systems  
– Length: 138 cm, height: 45 cm

# Mattresses

We consistently realize bock's aspiration to foster healthy rest. The lying surface is the key to comfort. Accordingly bock's mattresses fulfill special functions as a comfort component, adapted function intensifier or as a self-contained system for sleep. The focus is on health improvement in the choice of materials used in and around the mattresses. Special sizes are of course available on request.



- bock standard mattress RG 35/10 cm
  - PU-foam mattress, 90 x 200 x 10 cm
  - 1-side zip
    - (1 short side and 20 cm over the edge)
  - Cover Single Jersey (light blue), 90 g/m<sup>2</sup>, 100 % polyester
  - Upon request with PU-fitted sheet as a moisture protection
  - Alternative with PU-complete cover as a moisture protection



- bock standard mattress RG 40/12 cm
  - PU-foam mattress, 90 x 200 x 12 cm
  - 2-sided zip
    - (1 short side and 1 long side)
  - Cover made of elastic knitware, 180 g/m<sup>2</sup>, 58 % polyester, 42 % cotton
  - Washable 95 °C, autoclavable 105 °C
  - Upon request with PU-fitted cover as a moisture protection
  - Alternative with PU-complete cover as a moisture protection






- ▶ riposan-mattress RG 55 kg / m<sup>3</sup>
  - PU-foam mattress, approx. 14 cm high
  - 5-zone special anti-decubitus mattress
  - 2-sided zip
  - Suitable for allergic person
  - Dirt-resistant
  - Washable up to 95 °C



- ▶ Bariatric mattress XXL up to 200 kg
  - Foam mattress, approx. 20 cm high
  - Two outer cores made of cold foam and a highly resilient approx. 12 cm high special foam
  - Zipper on 4 sides,
  - High quality cover
  - Washable up to 60 °C



- ▶ Bariatric mattress with gel cover up to 300 kg
  - High lying comfort due to the special construction made of polyurethane, approx. 12 cm high
  - Zipper on two sides
  - High quality PU cover, blue
  - Washable up to 95 °C, suitable for dryer up to 50 °C

# New release 2016



## dino cot & adolescent care bed

New standards in care

### Specifications – dino cot & adolescent care bed

Total weight:	approx. 175 kg
Lying surface dimension:	<b>90 x 170 cm</b>
External dimension:	106.5 x 215 cm
Safe working load:	170 kg
Maximum person weight:	135 kg
Height adjustment:	16 – 86 cm
Lifter space:	> 15 cm
Side rail height :	74 cm
	Adptable to 109 cm and 144 cm

Total weight:	approx. 193 kg
Lying surface dimension:	<b>90 x 200 cm</b>
External dimension:	106.5 x 245 cm
Safe working load:	170 kg
Maximum person weight:	135 kg
Height adjustment:	16 – 86 cm
Lifter space:	> 15 cm
Side rail height :	74 cm
	Adptable to 109 cm and 144 cm

### Motors

Water-resistant to IPX4, high quality spiral cable with pull relief and anti-kink device.

### Features

Art.-No.	Lying surface dimension	External dimension
286	90 x 170 cm	106.5 x 215 cm
286	90 x 200 cm	106.5 x 245 cm

Subject to technical modification.



### 4-section

Electrical adjustment of the lying surface including the automatic triple knee break action.

### Basic equipment:

Bed frame 74 cm high  
Double folding doors with wooden pins and acrylic glass  
Electrical height adjustment  
Lying surface with removable aluminium bars  
High-stable 3-fold telescopic column system  
Solid wood in natural beech colour (clear varnished)  
Resistant to perspiration and spittle in accordance with DIN 71-3  
Safety stops at 22 cm and 34 cm

### Extras:

Height extension to 109 cm and 144 cm as retrofitting

### Accessories:

Lifting pole with hand grip (only for 90 x 200 cm)  
Interior bumpers  
Mattresses

### Safety guarantees:

- Primary fuse in the power plug
- Thermal release
- Locking function on the hand control
- Hand control with single faulty safety
- Special spiral mains cable
- Pull relief at power plug and anti-kink device
- Water-resistant to IPX4



Bed frame with view for ideal safety,  
optionally in the height of sides 74 cm, 109 cm or 144 cm

Maximum person weight 135 kg

92 cm (basic version )  
127 cm (high bed version)  
162 cm (high bed version)

Height adjustment  
16 - 86 cm

Doors can be folded away completely

Windows with acrylic glass

Quality lying surface with  
removable aluminium bars

Ø 100 mm castors with single  
brakes

External length of 215 cm for a lying surface length of 170 cm  
External length of 245 cm for a lying surface length of 200 cm

# New release 2016



## domiflex 2 wash with Push-and-Ready

Machine washable for the professional use

### Specifications – domiflex 2 wash

Total weight:	approx. 93 kg
Lying surface dimension:	90 x 200 cm
External dimension:	103 x 213 cm
Safe working load:	190 kg
Maximum person weight:	155 kg
Height adjustment:	39 – 80 cm
Lifter space:	> 15 cm
Side rail height :	39 cm
Special dimensions:	220 cm only width bed extension Special widths not available

### Motors

Water-resistant to IPX6, high quality spiral cable with pull relief and anti-kink device.

### Features

Art.-No.	Lying surface dimension	External dimension
272	90 x 200 cm	103 x 213 cm

Subject to technical modification.



### 4-section

Electrical adjustment of the lying surface including the automatic triple knee break action.



### Basic equipment:

Safety extra-low voltage system  
Adapted design for machine washing  
Lifting pole with hand grip  
Colour: F 37/009 Steamed Beech  
End panels with special infill  
Steel lying surface with removable bars  
Integrated steel side rails  
Push-and-Ready connection technique  
Transport and storage system  
Integrated Van-transport system

### Extras:


Bed extension for the length to 220 cm  
Comfort-sitting position / Reverse Trendelenburg

### Accessories:

Infusion holder  
Mattresses  
Mattress extension for the length to 220 cm  
Side rail attachment  
Splash water protection for the lifting motors  
Side rail bumpers  
Side rail bumpers for side rails with attachment  
And many more


### Safety guarantees:

- Primary fuse in the power plug
- Secondary fuse/polyswitch
- Thermal release
- Locking function on the hand control
- Hand control with single faulty safety
- Special spiral mains cable
- Pull relief at power plug and anti-kink device
- Water-resistant to IPX6
- Distortion-free, highly stable side rails



# Staying. Living. Careing.

In variable design for high demands




## We design rooms. Residential care

Order-No.: 890.01952 (German)  
Order-No.: 890.01977 (English)  
Order-No.: 890.01978 (French)  
Order-No.: 890.01979 (Italian)




## Chairs and tables

Order-No.: 890.01588 (German)




## Wardrobes and furniture

Order-No.: 890.01628 (German)



## .book freecare

Order -No.: 890.01997 (German)



## practico ultra low 9.5/80

Order-No.: 890.02142 (German)  
Order -No.: 890.02141 (English)

# www.bock.net

Visit us on the internet!

The latest news about our company and our products can be found on [www.bock.net](http://www.bock.net) and in the social network:



## **www.bock.net**

Besides our company and product presentation, you can find here the latest notes, appointments and many more information.

The configurator for our nursing care beds helps you to choose between the different variations, options and extras - experience how versatile bock is!



## **Facebook**

The fast overview of all news about our company and our products!



## **Youtube**

Besides our corporate film, you can find here interviews and assembly films in different languages.

Possible special colours for belluno and combiflex



► Black Forest  
light wMaple  
F04 / 054



► Bloomed Beech  
F37 / 023



► Bloomed Cherry  
planked  
F35 / 005



► Dark Walnut  
F26 / 019



► Albany Walnut dark  
F26 / 031

**Further information?**

Further information about the home care product range can be found easily on the internet. Or do you look for direct advice? You can reach us easily by phone or e-mail.

Our dedicated team is looking forward to help you!

**Perfectly suited for care.**

Hermann Bock GmbH  
Nickelstraße 12  
33415 Verl  
Germany

Phone: + 49 (0) 5246 9205 0  
Fax: + 49 (0) 5246 9205 25  
Internet: [www.bock.net](http://www.bock.net)  
E-Mail: [info@bock.net](mailto:info@bock.net)

® registered trade mark

Order no.: 890.01230

Issued: 04 | 2016

Subject to technical modification